OSYlink AS-CCS, Best.-Nr. 77 0870 00

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i Note! The OSYlink module and connected room terminals must be configured using the SystemOrganizer of the Flamenco system.

Caution! Electrostatic sensitive components

The printed circuit board includes electrostatic sensitive components. Therefore avoid touching.

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A2: Wall mounting

- Push out the three fixing clips 2 until the holes for the screws 4 appear.
- Screw the OSYlink 1 to the wall.

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
C System Status LEDs

OSYlink AS-CCS or OSYlink AS-VKS:


Flashing slowly (1s / 1s) = Connection to OSYnet is ok.
Steady light = No connection to OSYnet.

Flashing slowly (1s / 1s) = Speech line is operational.
Flashing quickly = Active speech connection or announcement.

Flashing slowly (1s / 1s) = CCS-Bus is ok.
Flashing quickly = Minimum one fault on the CCS bus

Steady light = OSYlink AS is not operational! It is waiting for configuration data from the OSY-ControlCenter

OSYlink AS-L200:

Flashing slowly (1s / 1s) = Connection to OSYnet ok.
Steady light = No connection to OSYnet.

Flashing slowly (1s / 1s) = Ward bus WCB0 is ok.
Flashing quickly = Minimum one fault on the ward bus WCB0

Steady light = OSYlink AS is not operational! It is waiting for configuration data from the OSY-ControlCenter

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OSYlink AS-CCS, order no. 77 0870 00

Interface for connecting the room terminals of a CCS ward (CCS bus) to the group bus (OSYnet) of a Flamenco nurse call system. Current consumption: 90 mA. For connecting a half-duplex speech system (CCS 1080 W) additionally a speech amplifier, order no. 74 8000 00, with the power supply unit, order no. 21 8000 00, is needed.

OSYlink AS-VKS, order no. 77 0871 00

Interface for connecting the room terminals of a VKS ward (CCS bus) to the group bus (OSYnet) of a Flamenco^{SE} call system. Current consumption: 90 mA. Additionally a speech amplifier, order no. 74 8000 00, with the power supply unit, order no. 21 8000 00, is needed.

OSYlink AS-L200, order no. 77 0872 00

Interface for connecting the room terminals of a L200 ward (ward bus WCB0) to the group bus (OSYnet) of a Flamenco nurse call system. Current consumption: 90 mA.

Note! The OSYlink module and connected room terminals must be configured using the SystemOrganizer of the Flamenco system.

Caution! Electrostatic sensitive components
The printed circuit board includes electrostatic sensitive components. Therefore avoid touching.

A Mounting

Optionally mounting on supporting rail or wall mounting in distribution boxes or technical rooms.

A1: Mounting on supporting rail (35 mm)

Click the OSYlink 1 onto the supporting rail 3.

A2: Wall mounting

1. Push out the three fixing clips 2 until the holes for the screws 4 appear.
2. Screw the OSYlink 1 to the wall.

B Setting the 10 pole coding switch

Setting the address (1 – 110)

For setting the address for the OSYlink on the group bus OSYnet use the code switches 1 to 7 with the values 1, 2, 4, 8, 16, 32, 64. Set the appropriate switch to ON.

Examples: Address 1 is set by using the code switch 1. Address 4 by the code switch 3. Address 24 is selected through addition: coding switch 4 and 5 (8 + 16 = 24).

Setting the connected system (according to figure B)

For setting the connected system use code switch 8 „W/G“ (W = simplex speech, G = duplex speech) and code switch 9 „1080/2000“.

Reset

To reset the OSYlink you have to set code switch 10 to ON for one second, then back to OFF. A reset has to be made, if a code switch setting has been changed. Otherwise the new setting will not be activated.

C System Status LEDs

OSYlink AS-CCS or OSYlink AS-VKS:

CCS L200 Audio OSYnet
green - green red
Flashing slowly (1s / 1s) = Connection to OSYnet is ok.
Steady light = No connection to OSYnet.
Flashing slowly (1s / 1s) = Speech line is operational.
Flashing quickly = Active speech connection or announcement.
Flashing slowly (1s / 1s) = CCS-Bus is ok.
Flashing quickly = Minimum one fault on the CCS bus
Steady light = OSYlink AS is not operational! It is waiting for configuration data from the OSY-ControlCenter

OSYlink AS-L200:

CCS L200 Audio OSYnet
- green - red
Flashing slowly (1s / 1s) = Connection to OSYnet ok.
Steady light = No connection to OSYnet.
Flashing slowly (1s / 1s) = Ward bus WCB0 is ok.
Flashing quickly = Minimum one fault on the ward bus WCB0
Steady light = OSYlink AS is not operational! It is waiting for configuration data from the OSY-ControlCenter

D Dismantling from supporting rail

Put a screw driver into the visible hole of the black fixing clip. Then pull the clip upwards, until the OSYlink 1 comes loose from the supporting rail 3.

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© Tunstall GmbH, Orkotten 66, D-48291 Telegate, www.tunstall.de
© Tunstall GmbH, Okoliten 66, D-48291 Telgte, www.tunstall.de