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Connections

1. Strip the connecting cables in the back box to a suitable length.
2. Remove the cable shield and drain wire up to the cable sheath. Strip the wires.
3. Connect the wires to the two connectors according to fig. **E**.

RAN-TEST (DIP switch)

Fig. **C**: Test if all room devices are ready to operate and if they are correctly connected to the room bus (RAN):

1. Set DIP switch “RAN-Test” to ON (left).
The LEDs of all connected room devices must flash. (Non-flashing devices are wrongly connected or defective. Replace defective devices. Correct any wrong connection.)
2. Compare the number of connected room devices with the RAN number configured with the ControlTerminal ConfigSet. (If the configured RAN number is not equal to the number of connected room devices, you have to correct the RAN number setting with the ControlTerminal ConfigSet.)
3. To end the test, set DIP switch “RAN-TEST” back to OFF (right).

Technical data

Power supply	24 V DC
Standby current consumption	42 mA
Current consumption per light section	30 mA
4-pole connector: - Connection capacity - Insulation strip length	0.20 – 2.50 mm ² 7 mm
8-pole connector: - Connection capacity - Insulation strip length	0.14 – 1.50 mm ² 7 mm
Dimensions (HxWxD)	190 x 150 x 40 mm
Field for the label strip (HxW)	70 x 92 mm
Degree of protection	IP 20
Ambient temperature	+5 °C – +40 °C
Relative humidity	0 % – 85 %

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
Mounting

At the point of delivery the ControlTerminal with doorplate is assembled and must be dismantled as follows, see fig. **A**:

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Wall mounting, see fig. **B**:

1. Connect the connection wires according to chapter “Connections” of this product leaflet to the connectors (screw clamps) on the mounting plate [2]. Thereby insert two insertion bridges as shown.
2. Unplug the four LED modules [6] and lay them aside.
3. Fit the mounting plate [2] to the back box [1] using the four back box screws [3]. For the correct mounting direction the connector for power supply must be at the top.
4. Place the casing [4] (incl. the pre-mounted PCB) onto the mounting plate [2]. Press very carefully, until plugs and sockets of the connectors are firmly connected.
5. Fit the casing [4] (incl. the pre-mounted PCB) with the two fixing screws [5] to the mounting plate [2].
6. Plug the four LED modules onto the PCB according to the fig. **D**.
7. Press the light dome with the insert module [7] onto the casing [4] until it locks in place.

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- 1 - * Back box
- 2 - ** Mounting plate with connectors
- 3 - * Four back box fixing screws
- 4 - Casing incl. pre-mounted PCB
- 5 - ** Two fixing screws
- 6 - Four LED modules

- 7 - Light dome with insert module
 - 8 - * Label strip
 - 9 - Protection cover
- Order separately:**
- * Back box (incl. four screws)
 - ** Mounting kit, order no. 77 0960 00

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0V — Einlegebrücke (00 0220 52)
 0V —

Power supply
 Cable: NYM 2x2x2.5 mm²

+24V — Insertion bridge (00 0220 52) from previous room or to next room
 +24V —

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Room bus (RAN) without speech
 Cable: IY(ST)Y 2x2x0.8 (max. 4 wires per connection point)

+24V-Si — red
 RAN — yellow
 0V — black
 0V — white
 to RAN user

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
* **NOTE!** The NF wires exist only at the group bus OSYnet, if ComTerminals (= terminals with speech) are connected to the same group bus. These wires are not connected to the ControlTerminal. Each of them has to be connected through using the push-wire connector (order no. 00 0222 88, wire diameter 0.6 – 0.8 mm).

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ControlTerminal with doorplate Flamenco, order no. 77 0551 00

Control and monitoring of all room functions in compliance with DIN VDE 0834. Direct connection to the group bus (OSYnet) and the room bus (RAN).

Optical signalling of all call types and staff presences as well as additional display for WC call. Doorplate as label field for room designation.

NOTE! Prior to mounting the ControlTerminal has to be configured using the ControlTerminal ConfigSet, order no. 77 0920 00. The software in the ConfigSet must have revision **1.04** or higher!

NOTE! The complete installation of the system is described in the Technical Manual.

CAUTION! The printed circuit board and the LED modules include electrostatic sensitive components. Avoid touching.

Connections

1. Strip the connecting cables in the back box to a suitable length.
2. Remove the cable shield and drain wire up to the cable sheath. Strip the wires.
3. Connect the wires to the two connectors according to fig. E.

RAN-TEST (DIP switch)

Fig. C: Test if all room devices are ready to operate and if they are correctly connected to the room bus (RAN):

1. Set DIP switch "RAN-Test" to ON (left).

The LEDs of all connected room devices must flash. (Non-flashing devices are wrongly connected or defective. Replace defective devices. Correct any wrong connection.)

2. Compare the number of connected room devices with the RAN number configured with the ControlTerminal ConfigSet. (If the configured RAN number is not equal to the number of connected room devices, you have to correct the RAN number setting using the ConfigSet.)
3. To end the test, set DIP switch "RAN-TEST" back to OFF (right).

Technical data

Power supply	24 V DC
Standby current consumption	42 mA
Current consumption per light section	30 mA
4-pole connector:	
- Connection capacity	0.20 – 2.50 mm ²
- Insulation strip length	7 mm
8-pole connector:	
- Connection capacity	0.14 – 1.50 mm ²
- Insulation strip length	7 mm
Dimensions (HxWxD)	190 x 150 x 40 mm
Field for the label strip (HxW)	70 x 92 mm
Degree of protection	IP 20
Ambient temperature	+5 °C – +40 °C
Relative humidity	0 % – 85 %

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Mounting


At the point of delivery the ControlTerminal with doorplate is assembled and must be dismantled as follows, see fig. A:

- Lightly compress the light dome [7] from top and bottom. Then pull off the light dome [7].


Wall mounting, see fig. B:

- Connect the connection wires according to chapter "Connections" of this product leaflet to the connectors (screw clamps) on the mounting plate [2]. Thereby insert two insertion bridges as shown.
- Unplug the four LED modules [6] and lay them aside.
- Fit the mounting plate [2] to the back box [1] using the four back box screws [3]. For the correct mounting direction the connector for power supply must be at the top.
- Place the casing [4] (incl. the pre-mounted PCB) onto the mounting plate [2]. Press very carefully, until plugs and sockets of the connectors are firmly connected.
- Fit the casing [4] (incl. the pre-mounted PCB) with the two fixing screws [5] to the mounting plate [2].
- Plug the four LED modules onto the PCB according to the fig. D.
- Press the light dome with the insert module [7] onto the casing [4] until it locks in place.

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
RESET: To reset the ControlTerminal you have to set DIP switch RESET to ON for one second, then back to OFF.

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Risk of short circuit! When setting the DIP switch do not contact any other electronic components.

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LED modules

- red: Calls
- green: Staff 1
- yellow: Staff 2
- white: WC call

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







- * Back box
- ** Mounting plate with connectors
- * Four back box fixing screws
- Casing incl. pre-mounted PCB
- ** Two fixing screws


- Four LED modules
- Light dome with insert module
- * Label strip
- Protection cover

Order separately:

- * Back box (incl. four screws)
- ** Mounting kit, order no. 77 0960 00

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*** NOTE!** The NF wires exist only at the group bus OSYnet, if ComTerminals (= terminals with speech) are connected to the same group bus. These wires are not connected to the ControlTerminal. Each of them has to be connected through using the push-wire connector (order no. 00 0222 88, wire diameter 0.6 – 0.8 mm).

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© Tunstall GmbH, Orkotten 66, 48291 Teigte, Germany, www.tunstall.de