**OSYlink-Universal, Best.-Nr. 77 0803 00**

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**OSYlink-Universal, order no. 77 0803 00**

Interface for connecting external systems and/or technical installations to the group bus (OSYnet).

Standby current consumption (all outputs OFF): 40 mA.

Max. current consumption (all outputs ON): 100 mA.

**Note!** The complete installation of the system is described in the technical manual.

**Caution!** The printed circuit board includes electrostatic sensitive components. Avoid touching.

**A Mounting**

Optionally mounting on supporting rail or wall mounting.

**A1: Mounting on supporting rail (35 mm)**

Click the OSYlink 1 onto the supporting rail 3.

**A2: Wall mounting**

- Push out the three fixing clips 2 until the holes for the screws 4 appear.
- Screw the OSYlink 1 to the wall.

**B Setting of address (1 - 110)**

For setting the address on the printed circuit board please use the code switches 1 to 7 with the values 1, 2, 4, 8, 16, 32, 64. Set the appropriate switch to ON.

Examples: Address 1 is set by using the code switch 1. Address 4 by the code switch 3.

Address 24 is selected through addition: coding switch 4 and 5 ( $8 + 16 = 24$ ).

Do not change code switches 8 and 9! Factory set: OFF.

**Reset**

To reset the OSYlink you have to set code switch 10 to ON for one second, then back to OFF. A reset has to be made, if a code switch setting has been changed. Otherwise the new setting will not be activated.


**C ESD protection cover (order no. 00 0276 54)**


After you have connected the cables (see reverse side), you have to put the ESD protection cover on the OSYlink. (ESD = Electro Static Discharge)

**D Dismantling from supporting rail**

Put a screw driver into the visible hole of the black fixing clip. Then pull the clip upwards, until the OSYlink 1 comes loose from the supporting rail 3.

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<b>IN 1</b>	Monitored input - Call	Normally open contact with parallel 4K7 resistor
<b>IN 2</b>	Monitored input - Emergency call	Normally open contact with parallel 4K7 resistor
<b>IN 3</b>	Monitored input- Alarm call	Normally open contact with parallel 4K7 resistor
<b>IN 4</b>	Monitored input - Call	Normally open contact with parallel 4K7 resistor
<b>Function for IN 1 - IN 4:</b> Contact closed = Call message Fresh call Contact opened = Off message Call cancelled Acknowledgement possible, fresh / answered		

<b>IN 5</b>	Input - Collective announcement (to all wards)	Normally open contact
<b>IN 6</b>	Input - Collective announcement (to all staff)	Normally open contact
<b>Function for IN 5 - IN 6:</b> Contact closed = Announcement started Contact opened = Announcement closed		

<b>IN 7</b>	Input - Raise call	Normally open contact
<b>IN 8</b>	Input - Cancel call	Normally open contact
<b>Function for IN 7:</b> Contact closed = Raise call Contact opened = no function		
<b>Function for IN 8:</b> Contact closed = Cancel call from IN7 Contact opened = no function		

<b>OUT 1</b>	Output can be configured. Factory setting: Collective output Calls	Solid state output 24V
<b>OUT 2</b>	Output can be configured. Factory setting: Collective output Emergency calls	Solid state output 24V
<b>OUT 3</b>	Output can be configured. Factory setting: Collective output Alarm calls	Solid state output 24V
<b>OUT 4</b>	Output can be configured	Solid state output 24V
<b>OUT 5</b>	Output can be configured	Solid state output 24V
<b>OUT 6</b>	Output can be configured	Solid state output 24V and/or location light feature for inputs IN 1, IN 2, IN 3, IN 4, IN 7.
<b>OUT 7</b>	Out. can be configured	Change-over contact potential-free, internal 24V can be used
<b>OUT 8</b>	Out. can be configured	Change-over contact potential-free, internal 24V can be used

The outputs are configured using the SystemOrganizer software.  
 In case of ward coupling the outputs signal across wards.

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### A2: Wall mounting

1. Push out the three fixing clips 2 until the holes for the screws 4 appear.
2. Screw the OSYlink 1 to the wall.

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Examples: Address 1 is set by using the code switch 1. Address 4 by the code switch 3. Address 24 is selected through addition, coding switch 4 and 5 ( $8 + 16 = 24$ ).

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

## D Dismantling from supporting rail

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<b>IN 1</b>	Monitored input - Call	Normally open contact with parallel 4K7 resistor
<b>IN 2</b>	Monitored input - Emergency call	Normally open contact with parallel 4K7 resistor
<b>IN 3</b>	Monitored input- Alarm call	Normally open contact with parallel 4K7 resistor
<b>IN 4</b>	Monitored input - Call	Normally open contact with parallel 4K7 resistor
<b>Function for IN 1 - IN 4:</b> Contact closed = Call message Fresh call Contact opened = Off message Call cancelled Acknowledgement possible, fresh / answered		
<b>IN 5</b>	Input - Collective announcement (to all wards)	Normally open contact
<b>IN 6</b>	Input - Collective announcement (to all staff)	Normally open contact
<b>Function for IN 5 - IN 6:</b> Contact closed = Announcement started Contact opened = Announcement closed		
<b>IN 7</b>	Input - Raise call	Normally open contact
<b>IN 8</b>	Input - Cancel call	Normally open contact
<b>Function for IN 7:</b> Contact closed = Raise call Contact opened = no function		
<b>Function for IN 8:</b> Contact closed = Cancel call from IN7 Contact opened = no function		
<b>OUT 1</b>	Output can be configured. Factory setting: Collective output Calls	Solid state output 24V
<b>OUT 2</b>	Output can be configured. Factory setting: Collective output Emergency calls	Solid state output 24V
<b>OUT 3</b>	Output can be configured. Factory setting: Collective output Alarm calls	Solid state output 24V
<b>OUT 4</b>	Output can be configured	Solid state output 24V
<b>OUT 5</b>	Output can be configured	Solid state output 24V
<b>OUT 6</b>	Output can be configured	Solid state output 24V and/or location light feature for inputs IN 1, IN 2, IN 3, IN 4, IN 7.
<b>OUT 7</b>	Out. can be configured Change-over contact potential-free, internal 24V can be used	
<b>OUT 8</b>	Out. can be configured Change-over contact potential-free, internal 24V can be used	

The outputs are configured using the SystemOrganizer software.  
In case of ward coupling the outputs signal across wards.

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© Tunstall GmbH, Okolten 66, D-48291 Telgte, www.tunstall.de