**Flurdisplay Alpha 16, Best.-Nr. 77 0150 00**

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**Corridor Display Alpha 16, order no. 77 0150 00**

Display for alphanumeric indication of calls and general system information.

- 16 digits, alphanumeric
- Longer texts are automatically displayed as ticker
- Buzzer for signalling forwarded calls
- Loudspeaker, e.g. for announcements or signal tones
- Supply current: 300 mA (average)

**Note!** The complete installation of the system is described in the technical manual.

**A Setting of address (1 - 110)**

For setting the address please use the code switches 1 to 7 with the values 1, 2, 4, 8, 16, 32, 64. Set the appropriate switch to ON.

Examples: Address 1 is set by using the code switch 1. Address 4 is set by using the code switch 3. Address 24 is selected through addition: coding switch 4 and 5 (8 + 16 = 24).


After the address is changed while system operation a **reset** has to be made. For the reset set code switch 10 for one second to ON.

**B Mounting**

The corridor display is designed for mounting in corridors or duty rooms.

The standard accessory is used for ceiling mounting.

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
**C Connection**

The corridor display is directly connected to the group bus OSYnet. Use the free wire ends of the connection cable for on site connection.


**D Functional test**

When the corridor display is connected, a functional test should be carried out.  
 ● To start the functional test set code switch 9 "Test" to ON.  
 The corridor display will run through the test steps described in the following table. Each step takes approx. 2 seconds.

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
Display shows:	Test ok, if...	What is being tested?
<b>Start Test</b>	<i>Functional test starts.</i>	
<b>Tunstall GmbH Display Check 0123456789abcdef</b>	The 3 texts are displayed one after the other.	Optical display.
<b>Buzzer test</b>	<i>Buzzer test starts.</i>	
<b>Buzzer: low</b>	Low signal tone sounds.	Buzzer setting "low".
<b>Buzzer: medium</b>	Medium signal tone sounds.	Buzzer setting "medium".
<b>Buzzer: loud</b>	Loud signal tone sounds.	Buzzer setting "loud".
<b>Buzzer OFF</b>	<i>Buzzer is switched off.</i>	
<b>Speech relay 1 ON</b>	Low relay click can be heard.	Speech channel 1 of system.
<b>Speech relay 2 ON</b>	Low relay click can be heard.	Speech channel 2 of system.
<b>Speech relays OFF</b>	Relay click can be heard.	Speech channels are closed.
 (displayed as ticker)	Display without "holes". Holes show defective LEDs.	LED function.
<b>To stop the test: DIP switch 9 -&gt; OFF</b>	<i>To stop the test set code switch 9 to OFF.</i>	
<b>End of Test</b>	<i>Test ended.</i>	
<b>00:00 (time)</b>	<i>If OSYnet is connected: System time is displayed. Otherwise internal clock of corridor display starts at 00:00 h.</i>	

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
**Red control LED flashes slowly**  
 (1 sec ON, 1 sec OFF), if connection to OSYnet is ok. Otherwise there is a fault.

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**D Functional test**


When the corridor display is connected, a functional test should be carried out.

- To start the functional test set code switch 9 "test" to ON. The corridor display will run through the test steps described in the following table. Each step takes approx. 2 seconds.

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
Display shows:	Test ok, if...	What is being tested?
<b>Start Test</b>	<i>Functional test starts.</i>	
<b>Tunstall GmbH Display Check 0123456789abcdef</b>	The 3 texts are displayed one after the other.	Optical display.
<b>Buzzer test</b>	<i>Buzzer test starts.</i>	
<b>Buzzer: low</b>	Low signal tone sounds.	Buzzer setting "low"
<b>Buzzer: medium</b>	Medium signal tone sounds.	Buzzer setting "medium"
<b>Buzzer: loud</b>	Loud signal tone sounds.	Buzzer setting "loud"
<b>Buzzer OFF</b>	<i>Buzzer is switched off.</i>	
<b>Speech relay 1 ON</b>	Low relay click can be heard.	Speech channel 1 of system
<b>Speech relay 2 ON</b>	Low relay click can be heard.	Speech channel 2 of system
<b>Speech relays OFF</b>	Relay click can be heard.	Speech channels are closed.
(displayed as ticker)	Display without "holes". Holes show defective LEDs.	LED function
<b>To stop the test: DIP switch 9 -&gt; OFF</b>	<i>To stop the test set code switch 9 to OFF.</i>	
<b>End of Test</b>	<i>Test ended.</i>	
<b>00:00 (time)</b>	<i>If OSYnet is connected: System time is displayed. Otherwise internal clock of corridor display starts at 00:00 h.</i>	

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
Red control LED flashes slowly (1 sec ON, 1 sec OFF), if the connection to OSYnet is ok. Otherwise there is a fault.

77 0150 00  
Afficheur de couloir Alpha 16  
Corridor Display Alpha 16

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© Tunstall GmbH, Okolten 66, D-48291 Telgte, www.tunstall.de
© Tunstall GmbH, Okoliten 66, D-48291 Telgte, www.tunstall.de