Flurdisplay Alpha 16, doppelseitig, Best.-Nr. 77 0160 00

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


A Mounting

The corridor display is designed for mounting in corridors or duty rooms. The standard accessory is used for ceiling mounting.

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B Connection

The corridor display is directly connected to the group bus OSYnet. Use the free wire ends of the connection cable for on site connection.

C Setting of address (1 - 110)

Behind the protection cover at the side of one display there is a 10 part code switch. For setting the address please use the code switches 1 to 7 with the values 1, 2, 4, 8, 16, 32, 64. Set the appropriate switch to ON.

Examples: Address 1 is set by using the code switch 1. Address 4 is set by using the code switch 3. Address 24 is selected through addition: coding switch 4 and 5 (8 + 16 = 24).


After the address is changed while system operation a **reset** has to be made. For the reset set code switch 10 for one second to ON.

D Functional test

After the corridor display has been connected, a functional test should be carried out.

- To start the functional test set code switch 9 "Test" to ON.

The corridor display will run through the test steps described in the following table. Each step takes approx. 2 seconds. You have to carry out the test twice, i.e. check each display separately.

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
Display shows:	Test ok, if...	What is being tested?
Start Test	<i>Functional test starts.</i>	
Tunstall GmbH Display Check 0123456789abcdef	The 3 texts are displayed one after the other.	Optical display.
Buzzer test	<i>Buzzer test starts.</i>	
Buzzer: low	Low signal tone sounds.	Buzzer setting "low".
Buzzer: medium	Medium signal tone sounds.	Buzzer setting "medium".
Buzzer: loud	Loud signal tone sounds.	Buzzer setting "loud".
Buzzer OFF	<i>Buzzer is switched off.</i>	
Speech relay 1 ON	Low relay click can be heard. ¹⁾	Speech channel 1 of system.
Speech relay 2 ON	Low relay click can be heard. ¹⁾	Speech channel 2 of system.
Speech relays OFF	Relay click can be heard. ¹⁾	Speech channels are closed.
 (displayed as ticker)	Display without "holes". Holes show defective LEDs.	LED function.
To stop the test: DIP switch 9 -> OFF	<i>To stop the test set code switch 9 to OFF.</i>	
End of Test	<i>Test ended.</i>	
00:00 (time)	<i>If OSYnet is connected: System time is displayed. Otherwise internal clock of corridor display starts at 00:00 h.</i>	

¹⁾ The relays for the speech channels are only in the display with the coding switch.

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
Red control LED flashes slowly (1 sec ON, 1 sec OFF), if connection to OSYnet is ok. Otherwise there is a fault.

77 0160 00
Flurdisplay Alpha 16, doppelseitig
Corridor Display Alpha 16, double-sided

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Installation Instructions

GB

Corridor Display Alpha 16, double-sided, order no. 77 0160 00

Double-sided display for indication of calls and general system information.

- 16 digits, alphanumeric
- Longer texts are automatically displayed as ticker
- Tone generator for signalling forwarded calls
- Loudspeaker, e.g. for announcements or signal tones
- Supply current: 600 mA (average)




Note! The complete installation of the system is described in the technical manual.

A Mounting

The corridor display is designed for mounting in corridors or duty rooms.

The standard accessory is used for ceiling mounting.

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

Display shows:	Test ok, if...	What is being tested?
Start Test	<i>Functional test starts.</i>	
Tunstall GmbH Display Check 0123456789abcdef	The 3 texts are displayed one after the other.	Optical display.
Buzzer test	<i>Buzzer test starts.</i>	
Buzzer: low	Low signal tone sounds.	Buzzer setting "low"
Buzzer: medium	Medium signal tone sounds.	Buzzer setting "medium"
Buzzer: loud	Loud signal tone sounds.	Buzzer setting "loud"
Buzzer OFF	<i>Buzzer is switched off.</i>	
Speech relay 1 ON	Low relay click can be heard. ¹⁾	Speech channel 1 of system
Speech relay 2 ON	Low relay click can be heard. ¹⁾	Speech channel 2 of system
Speech relays OFF	Relay click can be heard. ¹⁾	Speech channels are closed.
 (displayed as ticker)	Display without "holes". Holes show defective LEDs.	LED function
To stop the test: DIP switch 9 -> OFF	<i>To stop the test set code switch 9 to OFF.</i>	
End of Test	<i>Test ended.</i>	
00:00 (time)	<i>If OSYnet is connected: System time is displayed. Otherwise internal clock of corridor display starts at 00:00 h.</i>	

¹⁾ The relays for the speech channels are only in one display with the coding switch.

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© Tunstall GmbH, Okoliten 66, D-48291 Telgte, www.tunstall.de