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
Corridor display, 16-digit, double-sided

The corridor display is intended for connection to the ward bus of a CONCENTO^{CARE} or CONCENTO^{PLUS} nurse call system. The power is supplied not via the ward bus but via a separate stub line from a power supply unit of the ward. The corridor display is 16-digit and displays calls with call type and call location.

NOTE! The complete installation of the system is described in the Technical Manual.

- [1] Flurdisplay
- [2] Anschlussleitung (ca. 3 m)
- [3] Lichtsensor
- [4] Dübel (Ø 8 mm)
- [5] Stockschraube
- [6] Deckenbefestiger

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- [7] Fitting
- [8] Screw cap
- [9] Wire rope (approx. 500 mm)
- [10] Head element of the wire rope holder
- [11] Wire rope holder

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Mounting

Ceiling suspension, not more than 2.5 m distance to junction box.

Attention! Take into account the light sensor position when aligning the corridor display.

1. Drill two holes in the ceiling (distance: 525 mm). Insert the dowels [4].
2. Each wire rope: Bolt together the wire rope [9] with fitting [7], ceiling fastener [6], screw cap [8] and hanger bolt [5] according to fig. C.
3. Each wire rope: Slide the wire rope [9] into the wire rope holder [11] up to the desired position. To lock the wire rope [9], pull it. A gripping mechanism is activated. Fig. D.

Note! To release an inadvertently activated locking, press the head element [10] of the wire rope holder [11]. The gripping mechanism will be released.

4. Screw the two hanger bolts [5] of the pre-mounted corridor display into the dowels [4] in the ceiling.

Optional accessories: Extension set ceiling suspension, order no. 19 0780 05, to extend the ceiling suspension by 50 cm.

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E Connection

The corridor display is looped into the ward bus with the pre-assembled connection cable and connected to a power supply unit of the ward with a NYM stub line.

1. The ends of the connection cable of the corridor display are tinned. Strip the other connecting cables in the junction box to a suitable length. Remove the cable shield and drain wire up to the cable sheath. Strip the wires.
2. Connect the wires to a connection box terminal or luster terminal according to the connecting diagram.

Note! The connecting cable of the corridor display can be shortened if necessary.

Bus termination

If the corridor display is the first or the last ward bus user, you have to set a 2.7 kOhm resistor between IA and IB at the connection point.


F Setting the address

In the factory settings the corridor display has the address: 9980. When two or more corridor displays are connected to a single ward bus, the second and further displays must be set to other addresses. Available addresses are from 9980 to 9989.

To set up the address, press the deep-set button (under the protection cover) at the top of the corridor display using a ballpoint pen. Press the button according to fig. F.

10 s = 10 seconds, 3 s = 3 seconds, 1x = Short-press once.

Timeout: If the button is not pressed for longer than 1 minute, the procedure is terminated automatically without saving.

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Technical data

Power supply	24 V DC
Current consumption for call	350 mA
Dimensions (HxWxD)	145 x 785 x 55 mm
Weight	Approx. 2500 g
Housing material	Lacquered aluminium
Front glass material	Acrylic glass
Degree of protection	IP 20
Ambient temperature	+5°C – +40°C
Relative humidity	0 % – 85 %