Anwesenheitstaster-Einsatz

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Presence switch insert

Intended for connection to the room bus of a CONCENTO^{CARE} nurse call system.


NOTE! The complete installation of the system is described in the Technical Manual.

A Product description

A green presence button serves for switching the presence of the staff on and off.

An LED in the presence switch lights up when the presence is switched on (reminder light).

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- [1] *Back box, 1-gang
(solid wall: 17 0100 00,
partition wall: 17 5100 00)
[2] Mounting plate
[3] *Back box screws

- [4] *Cover frame, 1-gang
[5] Terminal block
(plug-in terminals)
[6] Presence switch insert
[7] Presence button (green)

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* Not included in the scope of delivery, please order separately.

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B Mounting

Wall mounting on one-gang back box.

1. Bolt the mounting plate [2] to the back box [1]. The four square recesses in the mounting plate are located at the top and bottom.
2. Prepare the wires and feed them through the mounting plate [2] (see section "Connection").
3. Connect the switch insert [6] according to section "Connection".
4. Hold the cover frame [4] in front of the mounting plate [2] on the wall.
5. Plug the presence switch insert [6] onto the mounting plate so that the latching hook of the switch insert fits into the square recesses of the mounting plate.

C Dismantling


1. Pull the cover frame [4] together with the switch insert [6] off the mounting plate [2].
2. For disconnecting the wires, refer to fig. E.

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
EN - Installation Instructions

D Connection

1. Pull the connection cable out through the mounting plate and strip it to a suitable length. Remove the cable shield and drain wire up to the cable sheath. Strip the wires.
2. Connect the wires according to the connection diagram.

Note: The terminal block [5] can be pulled off temporarily.

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Technical data

Nominal voltage	24 V=
Standby current consumption	8 mA
Room bus connection	Cable type: J-Y(St)Y 2x2x0.8
Skinning length	6 mm
Installation depth	24 mm
Degree of protection	IP 22
Ambient temperature	+5°C – +40°C
Relative humidity	0 % – 95 %