## Flurdisplay, 12-stellig

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
## Corridor display, 12-digit

The corridor display is intended for connection to the ward bus of a CONCENTO<sup>CARE</sup> nurse call system. The power is supplied not via the ward bus but via a separate stub line from a power supply unit of the ward.

The corridor display is 12-digit and displays calls with call type and call location.

**NOTE!** The complete installation of the system is described in the Technical Manual.

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- [1] Corridor display  
[2] Connection cable (approx. 3 m)

- [3] Light sensor

### A Wall mounting

1. Hang the corridor display on the wall with the aid of the keyhole-shaped openings on the rear (distance to the junction box < 2.5 m).
2. Connect the wires of the pre-assembled connection cable according to section "Connection".

### B Dismantling

1. Open the junction box and disconnect the connections.
2. Remove the corridor display from the suspension.

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**C Connection**

The corridor display is looped into the ward bus with the pre-assembled connection cable and connected to a power supply unit of the ward with a stub line.

1. The ends of the connecting cable of the corridor display are tinned. Strip the other connecting cables in the junction box to a suitable length. Remove the cable shield and drain wire up to the cable sheath. Strip the wires.
2. Connect the wires to a connection box terminal or luster terminal according to the connecting diagram.

**Note!** The connecting cable of the corridor display can be shortened if necessary.

**D Setting the address**

In the factory settings the corridor display has the address: 9980. When two or more corridor displays are connected to a single ward bus, the second and further displays must be set to other addresses. Available addresses are from 9980 to 9989.

To set up the address, press the deep-set button (under the protection cover) at the top of the corridor display using a ballpoint pen. Press the button according to fig. D.

10 s = 10 seconds, 3 s = 3 seconds, 1x = Short-press once.

**Timeout:** If the button is not pressed for longer than 1 minute, the procedure is terminated automatically without saving.

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**Technical data**

Power supply	24 V DC
Current consumption for call:	300 mA
Dimensions (HxWxD)	125 x 604 x 55 mm
Weight	1600 g
Housing material	Lacquered aluminium
Front glass material	Acrylic glass
Degree of protection	IP 20
Ambient temperature	+5°C – +40°C
Relative humidity	0 % – 85 %